



Voluson Performance Series








PROBE GUIDE

Voluson™

A Healthier Future for Women

gehealthcare.com



	Description	Applications	Bandwidth	FOV
CURVED - 2D				
 4C-RS H4000SR	Wideband Convex Probe	Abdomen, Obstetrics, Gynecology	2 – 5 MHz	81°
CURVED - REAL-TIME 4D				
 RAB2-6-RS H48681WR	Wideband Convex Ultra-light Volume Probe	Abdomen, Obstetrics, Gynecology	2 – 5 MHz	90° V 90° x 85°
ENDOCAVITY - 2D				
 IC9-RS H48691PJ	Wideband Microconvex Endocavitary Probe	Obstetrics, Gynecology, Transrectal	2.9 – 9.7 MHz	185°
 IC9B-RS H48062AK	Wideband Microconvex Endocavitary Probe	Obstetrics, Gynecology, Transrectal	2.9 – 9.7 MHz	185°
ENDOCAVITY - REAL-TIME 4D				
 RIC5-9A-RS H48701EJ	Wideband Microconvex Endocavitary Volume Probe	Obstetrics, Gynecology, Transrectal	3.8 - 9.3 MHz	184° V 184° x 120°
LINEAR - 2D				
 12L-RS H40402LY	Wideband Linear Probe	Small Parts, Peripheral Vascular, Pediatrics, Musculoskeletal, Breast	4 – 12 MHz	38 mm
PHASED ARRAY - 2D				
 3SC-RS H45041DL	Wideband Phased Array Probe	Abdominal, Obstetrics, Cardiology, Pediatrics, Cephalic	1 – 4 MHz	90°

For probe care and cleaning information,
visit www.gehealthcare.com/transducers.

© GE, 2021

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram, and Voluson are trademarks of GE. GE Healthcare, a division of GE. GE Medical Systems, Inc., doing business as GE Healthcare.

October 2021
JB17424XX

